

Able2

Rubber Ring

PR20616 / PR20618

60618 60618



Please inspect all parts carefully
for cracks/breaks prior to use.



Do not over-inflate

Distributed By:

Able2 Europe B.V.
Koperslagerij 15
4762 AR Zevenbergen
The Netherlands
www.able2.nl

Able2 UK Ltd.
www.able2.eu

Able2 Iberia S.L.
www.able2.es

USER INSTRUCTIONS


- Features special valve for easy inflation
- Twist open top of valve
- Inflate by mouth or bicycle-style pump to maximum height of 3 inches (7.6 cm)
- Tighten valve to close

CARE INSTRUCTIONS

- Wipe with damp cloth



Maximum
User Wt:



Valentine International Ltd.
8th Floor, No. 149, Sec 2, Ta Tung Road
221 Hsichih District, New Taipei City
Taiwan R.O.C.



Y. Sung Handelsvertretung
Duesselthaler Str. 24
40211 Duesseldorf Germany

Made in China



Rev. 071520